

UNIVERSITI  
KEBANGSAAN  
MALAYSIA

---

*The National University  
of Malaysia*

User's Manual


For

Labelme


## Introduction:


Despite training the model extensively on the given dataset, there are weaknesses to this model and therefore it may predict the wrong cell type. In such cases you may enter your own expert opinion using a built-in functionality of the app. This guide aims to get users familiar with the step-by-step process of labelling images.


## Scenario:




Solution:







labelme - C:/Users/同/Desktop/HA200041\_JPG.rf.34c9e97ddee69612f22c16cc2cc7900c.jpg



Navigate to:  
Edit→CreateRectangle  
And select CreateRectangle.



Flags

Label List


Polygon Labels

File List

Search Filename

Start drawing rectangles





Then go to the cell you want to make corrections to, and then draw a


**TIGHT FITTING** box around it.


To draw a box:

Click once on starting point A.

Click once on ending point B.

You now have a box with diagonals A and B.





After completing the box, a popup will appear. Enter the correct name of the cell in that space.

**Please make sure to keep the cell name consistent. E.g. neutrophil and Neutrophil are two very different things to the model, and labelling the data with incorrect or inconsistent spellings/capitalization will cause confusion.**

These are pre-existing cell types that you can use:

Lymphoblast

Neutrophil

Promyelocytes

Dyplastic\_Cell

Lymphocytes

myelocytes


Band\_Neutrophils

nucleated\_RBC


Giant\_Platelet

monocytes

alternatively, you can add new cell types not present in this list. Simply type in the new cell name and it will be added automatically to the cell types list.







After which you will select the folder to save your data into.

END



